

AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/20

AFPW00620





ADDITIONAL INFORMATION

Weight 0.77 kg

Dimensions $15 \times 3.5 \times 10 \text{ cm}$

Carton Qty 36

Pricing Unit Box/20

Size 6cm x 7cm

Min Order Qty 1

ARTG ID 159198 **GMDN Code** 44990

Medical Device Class 1s

UNSPSC 42311527

Tax Schedule Non-Taxable

Item Length 7
Item Width 6

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

AeroFilm Plus Waterproof Island Dressing provides transparent, waterproof protection with a soft absorbent island pad for comfortable post-operative and delicate wound coverage. The clear film acts as a bacterial barrier while allowing moisture vapour to escape, helping the site stay protected in daily life and letting users visually check the wound without removing the dressing.

The absorbent SupaSorb pad draws light exudate away from the skin and the low-adherent contact layer minimises trauma at change time. This combination supports a cleaner wound environment and gentler, low-pain removal—ideal for cuts, grazes, sutures and small procedural sites.

SupaBond hypoallergenic adhesive delivers reliable hold over body contours and joints, resisting edge lift from friction, sweat and movement. The film is ultra-thin for a discreet, conformable fit that moves with the body, maintaining a secure seal without bulk.

A unique frame delivery system simplifies accurate placement. The frame helps apply the dressing crease-free and ensures firm edge contact for a consistent seal in clinical settings and at home.



AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/20

AFPW00620

- Waterproof transparency: clear film shields the site while enabling visual checks.
- **Comfort absorbency:** SupaSorb pad manages exudate to help reduce maceration.
- **Hypoallergenic hold:** SupaBond adhesive provides strong, skin-friendly adhesion.
- Low-adherent contact: minimises sticking for gentler dressing changes.
- Easy, precise placement: frame delivery supports crease-free application.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

CONTENTS LISTS

Code	Description	Qty
AFPW006	AEROFILM PLUS Waterproof Island Dressing 6 x 7cm	1

FILM DRESSING PAD