



ADDITIONAL INFORMATION

Weight	0.77 kg
Dimensions	15 × 3.5 × 10 cm
Carton Qty	36
Pricing Unit	Box/20
Size	6cm x 7cm
Min Order Qty	1
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC Code	42311527
Tax Schedule	Non-Taxable
Item Length	7
Item Width	6
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AeroFilm Plus Waterproof Island Dressing pairs a transparent, waterproof film with a soft, absorbent island pad to protect delicate and post-operative sites through everyday activity. The clear polyurethane-style film forms a barrier to water and external contaminants while allowing moisture vapour to pass, helping the skin feel fresher under the dressing. Because the site remains visible through the film, quick visual checks can be made without disturbing the seal. Each dressing is sterile and single-use, making it suitable for clinical care, workplace first aid stations, and active households.

At the centre, the SupaSorb absorbent pad manages light exudate to help reduce maceration. A low-adherent contact layer minimises sticking to fragile tissue so removal is gentler and less disruptive. This balanced design supports a cleaner environment around cuts, grazes, sutures and small procedural sites, while the low profile sits comfortably under clothing and gloves.

SupaBond hypoallergenic adhesive provides dependable hold over contours and joints. It resists edge lift from friction, perspiration and motion to help the dressing stay in place during showering, work shifts and daily routines. The ultra-thin, conformable film moves smoothly with the body for a discreet finish that supports confidence during wear. The dressing is latex free and suitable for sensitive



AEROFILM PLUS WATERPROOF ISLAND DRESSING

6 X 7CM BOX/20

AFPW00620

skin.

An integrated frame delivery system guides accurate placement. The frame helps apply the dressing crease-free and ensures firm edge contact for a consistent, reliable seal. This simple, intuitive application is valued in clinics for speed and repeatability and at home for ease of use and consistent results. From minor procedures to routine first aid, AeroFilm Plus delivers transparent, waterproof protection with comfortable absorbency.

- **Waterproof transparency:** clear film shields the site and enables visual checks.
- **Comfort absorbency:** SupaSorb pad manages exudate to help reduce maceration.
- **Hypoallergenic strong hold:** SupaBond adhesive provides reliable adhesion on contours.
- **Low-adherent contact:** minimises sticking for gentler, low-pain removal.
- **Easy, precise placement:** frame delivery supports smooth, crease-free application.

The **6 × 7 cm size** offers generous coverage for small post-operative sites and sutures while remaining comfortable on contoured areas. It balances a broad protective film border with a soft, absorbent island pad for everyday movement.

Packaged in a **quantity of 20** to streamline restocking for kits and clinics, providing clear waterproof coverage, a comfortable absorbent island pad and hypoallergenic adhesion for dependable daily wear.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Film Dressing Pad