



ADDITIONAL INFORMATION

Weight	0.17 kg
Dimensions	14.7 × 8 × 10 cm
Carton Qty	8
Pricing Unit	Box/50
Size	6cm x 7cm
Min Order Qty	1
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311527
Tax Schedule	Non-Taxable
Item Length	7
Item Width	6
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36

PRODUCT DESCRIPTION

This post-operative dressing combines the advantages of AEROFILM with a SupaSorb™ absorbent and low adherent dressing pad to provide a superior dressing for post operative and delicate wound sites.



AEROFILM PLUS WATERPROOF ISLAND DRESSING

6 X 7CM BOX/50

AFPW006

AeroFilm
Ultra-thin, Waterproof & Transparent Bacterial Barrier

Supabond™
Hypo-allergenic adhesive
exclusively engineered,
for low irritation and high adhesion

SuperSorb™ Wound
Exudate Pad
Super Absorption,
Reducing Skin Maceration

Non-stick Barrier
Protects Wound with a Non-Adherent Barrier,
Minimising Trauma on Removal

Dressing Frame
Ensures Easy, Effective Application

SupaBond™ Hypo-allergenic Adhesive

- Exclusively engineered for unrivalled adhesion, this hypoallergenic adhesive has a low irritation risk while ensuring the dressing stays in place over body contours.

Transparent AeroFilm

- The transparent AeroFilm allows the user to monitor the level of fluid and check for clinical signs of infection without removing the dressing. This allows for maximum usage of the dressing before removal.

Simple Application

- AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over the wound site.

Contraindications / Precautions



AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/50

AFPW006

- AeroFilm should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

CONTENTS LISTS

Code	Description	Qty
RES150230	Clear Resealable Bag 15 x 23cm (9" x 6")	1

FILM DRESSING PAD