

## AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/3

**AFPW0063** 





### **ADDITIONAL INFORMATION**

Weight 0.023 kg

**Dimensions**  $1.5 \times 9.5 \times 15 \text{ cm}$ 

Carton Qty 96

Pricing Unit Box/3

**Size** 6cm x 7cm

Min Order Qty 1

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311527

**Tax Schedule** Non-Taxable

Item Length7Item Width6

**Item Measurement Unit** cm

Is Serialized Item No
Is Dangerous Good No

Shelf Life (Months) 36

#### **PRODUCT DESCRIPTION**

This post-operative dressing combines the advantages of AEROFILM with a SupaSorb<sup> $\mathbf{m}$ </sup> absorbent and low adherent dressing pad to provide a superior dressing for post operative and delicate wound sites.

AeroFilm



## AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/3

**AFPW0063** 

Ultra-thin, Waterproof & Transparent Bacterial Barrier

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

SuperSorb™ Wound Exudate Pad Super Absorption, Reducing Skin Maceration

Non-stick Barrier Protects Wound with a Non-Adherent Barrier, Minimising Trauma on Removal

**Dressing Frame** 

Ensures Easy, Effective Application

## SupaBond™ Hypo-allergenic Adhesive

• Exclusively engineered for unrivalled adhesion, this hypoallergenic adhesive has a low irritation risk while ensuring the dressing stays in place over body contours.

## **Transparent AeroFilm**

• The transparent AeroFilm allows the user to monitor the level of fluid and check for clinical signs of infection without removing the dressing. This allows for maximum usage of the dressing before removal.

### **Simple Application**

AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the
dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over the wound
site.



# AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/3

**AFPW0063** 

### **Contraindications / Precautions**

- AeroFilm should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.

Film Dressing Pad