



## ADDITIONAL INFORMATION

<b>Weight</b>	0.104 kg
<b>Dimensions</b>	16 × 5 × 16 cm
<b>Carton Qty</b>	20
<b>Pricing Unit</b>	Box/10
<b>Size</b>	9 x 10cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	407191
<b>GMDN Code</b>	58301
<b>Medical Device Class</b>	Class 2a
<b>Tax Schedule</b>	Non-Taxable
<b>Item Length</b>	10
<b>Item Width</b>	9
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



## PRODUCT DESCRIPTION

AeroFilm Plus Sensitive Silicone Postoperative Dressing is designed to support comfortable care of acute wounds while being kind to fragile and sensitive skin. Each dressing pairs a soft absorbent pad with a transparent film that seals securely around the site, helping to protect the area from external moisture and everyday contact. The silicone adhesive layer adheres gently to skin, allowing the dressing to be lifted and repositioned with minimal discomfort when required.

The absorbent pad manages wound exudate to help prevent pooling and excess moisture around the site. This balance supports a moist environment without leaving the surrounding skin overly wet, which can contribute to maceration. The dressing is suitable for a wide range of acute wounds such as postoperative incisions, cuts and abrasions, making it a versatile choice for hospital wards, day

surgery, outpatient clinics and home recovery.

The thin, flexible film conforms closely to body contours, including areas that bend and move frequently, such as elbows and knees. This flexibility helps the dressing stay comfortably in place during normal daily activity. The waterproof outer layer provides an effective barrier against water from everyday routines, helping to maintain a clean environment around the wound between dressing changes.

The clear film keeps the dressing discreet on the skin and allows the surrounding area to remain visible so that changes in the condition of the skin can be monitored without fully removing the product. Combined with the cushioning pad, this construction offers both protection and comfort for patients recovering from surgery or managing short-term wounds in clinical or home settings.

- **Gentle silicone adhesion:** designed for fragile and sensitive skin with comfortable removal.
- **Effective exudate management:** absorbent pad helps control wound fluid for an optimal environment.
- **Waterproof protective film:** shields the site from external moisture and everyday contact.
- **Highly flexible design:** conforms to body contours for comfortable wear on moving areas.
- **Discreet, transparent appearance:** low-profile film keeps the dressing unobtrusive on skin.

The **Size 9 x 10 cm** option provides generous coverage for medium acute wounds, offering a balanced area of absorbent pad and adhesive border that sits comfortably on limbs or torso while remaining low profile.

Supplied in a **box of 10 single-use dressings**, this option offers a practical balance between availability and storage space for clinical areas and home recovery. The 9 x 10 cm format provides consistent coverage for medium acute wounds, helping streamline dressing changes without overstocking.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Silicone Postoperative Dressing 9 x 10cm