



ADDITIONAL INFORMATION

Weight	0.14 kg
Dimensions	19.7 × 3.5 × 14 cm
Carton Qty	8
Pricing Unit	Box/20
Size	10cm x 12cm
Min Order Qty	1
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC Code	42311527
Tax Schedule	Non-Taxable
Item Length	12
Item Width	10
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36



PRODUCT DESCRIPTION

AeroFilm is a sterile, waterproof, ultra-thin transparent film dressing engineered for dependable protection on clean skin, superficial wounds and I.V. devices. The breathable film forms an effective barrier to bacteria and viruses while permitting moisture vapour exchange to help maintain a cleaner environment around the site. Its flexible construction conforms smoothly to contours and sits discreetly under clothing, supporting comfortable wear throughout daily activity in clinical, workplace and at-home settings.



AEROFILM WATERPROOF FILM DRESSING 10 X 12CM BOX/20

AFW010

SupaBond hypoallergenic adhesive is designed for strong, consistent hold that resists edge lift from movement and moisture. The latex-free formulation supports sensitive skin and removes cleanly when it is time to change. Because the film is crystal clear, the application site can be observed continuously without disturbing the area, helping reduce unnecessary changeovers and supporting good practice monitoring for I.V. insertion sites and minor wound coverage.

An integrated frame applicator simplifies handling and placement for crease-free results with minimal contact to the adhesive surface, even on curved or awkward body contours. The low-profile, ultra-thin film flexes with the body to reduce catching on clothing and helps provide a smooth, low-friction outer surface. Each dressing is sterile and single-use to support hygiene requirements. AeroFilm provides reliable coverage for minor burns, blisters and grazes, and is equally suited to protecting post-procedure areas where a waterproof, breathable barrier is desired.

- **Waterproof barrier:** protects against water and body fluids while remaining breathable for comfort.
- **Secure, skin-friendly hold:** SupaBond hypoallergenic adhesive resists edge lift and is latex-free.
- **Clear visibility:** transparent film enables site checks without removal.
- **Easy placement:** integrated frame applicator aids smooth, crease-free application.
- **Conformable comfort:** ultra-thin film moves with the body and wears discreetly under clothing.

The **10 × 12 cm size** provides generous coverage for post-procedure areas and I.V. sites, offering a clear, waterproof, breathable barrier that conforms smoothly to curves and joints.

Supplied in a **box of 20 dressings** to support convenient replenishment for kits and clinical trolleys while reducing mid-week reorder frequency for busy workplaces.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Waterproof Film Dressing