

## AEROFILM WATERPROOF FILM DRESSING 10 X 12CM BOX/20

**AFW010** 





## **ADDITIONAL INFORMATION**

Weight 0.14 kg

**Dimensions**  $19.7 \times 3.5 \times 14 \text{ cm}$ 

Carton Qty 8

Pricing Unit Box/20

Size 10cm x 12cm

Min Order Qty 1

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311527

**Tax Schedule** Non-Taxable

Item Length12Item Width10

**Item Measurement Unit** cm

**Is Medical Device** Yes

**Is Serialized Item** No

**Is Dangerous Good** No

Shelf Life (Months) 36

## **PRODUCT DESCRIPTION**

AeroFilm is a sterile, waterproof, ultra-thin transparent film dressing designed for dependable protection on clean skin, superficial wounds, and I.V. devices. The breathable film acts as a barrier to bacteria and viruses while allowing moisture vapour exchange, helping maintain a cleaner environment around the application site. Its conformable, flexible construction moves with the body and sits discreetly under clothing, making it suitable for everyday use in clinical, workplace, and at-home settings.

SupaBond hypoallergenic adhesive is engineered for strong, consistent hold that resists edge lift during movement. The latex-free formulation supports sensitive skin and removes cleanly when it's time to change. The crystal-clear film enables continuous site observation without disturbing the area, supporting good practice and reducing unnecessary changeovers.

An integrated frame applicator simplifies placement for crease-free results and minimal handling of the adhesive surface, even on curved or awkward body contours. Each dressing is sterile and single-use. AeroFilm provides reliable coverage for minor burns, blisters



## AEROFILM WATERPROOF FILM DRESSING 10 X 12CM BOX/20

**AFW010** 

and grazes, and is also well suited to protecting I.V. insertion sites and post-procedure areas where a waterproof, breathable barrier is required.

- Waterproof barrier: protects against water and body fluids while remaining breathable.
- Secure, skin-friendly hold: SupaBond hypoallergenic adhesive resists edge lift and is latex-free.
- Easy, accurate placement: frame applicator helps achieve a smooth, crease-free seal.
- Clear visibility: transparent film allows site monitoring without removal.
- Sterile single-use: suitable for first aid kits, clinics and workplace cabinets.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Waterproof Film Dressing