

## AEROFILM WATERPROOF FILM DRESSING 6 X 7CM

AFW0063





### **ADDITIONAL INFORMATION**

Weight 0.272 kg

**Dimensions**  $19 \times 9.5 \times 15 \text{ cm}$ 

Carton Qty 96

Pricing Unit Box/3

Size 6cm x 7cm

Min Order Qty 1

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311527

**Tax Schedule** Non-Taxable

Item Length7Item Width6

**Item Measurement Unit** cm

Is Medical Device Yes
Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 36

### **PRODUCT DESCRIPTION**

AeroFilm is a sterile, waterproof, ultra-thin, transparent film dressing. It is designed to be applied over I.V. devices, clean skin or wounds and provides a breathable anti-bacterial and viral barrier to outside contaminants.



# AEROFILM WATERPROOF FILM DRESSING 6 X 7CM BOX/3

**AFW0063** 

| AeroFilm<br>Ultra-thin, Waterproof & Transparent Bacterial Barrier   |
|--|
| Supabond™<br>Hypo-allergenic adhesive<br>exclusively engineered,<br>for low irritation and high adhesion   |
| Dressing Frame<br>Ensures Easy, Effective Application  |
| SupaBond adhesive  |
| <ul> <li>Exclusively engineered for unrivalled adhesion, ensures that AeroFilm can be left in place for up to 7 days. This reduces interference with the wound, reducing nursing time and aiding faster recovery.</li> </ul>   |
| Waterproof Film  |
| <ul> <li>AeroFilm is impermeable to water and body fluids, protecting the application underneath and allowing the patient to<br/>shower and bathe without redressing. It also aids in the prevention of bacterial contamination, and reduces risks of skir<br/>damage through interference.</li> </ul> |
| Simple application   |



### AEROFILM WATERPROOF FILM DRESSING 6 X 7CM BOX/3


AFW0063

|   | AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over the wound site. |
|---|--|
| 0 | Awkward Areas  AeroFilm can be easily applied to the most awkward body contours, due to its expandability, while SupaBond ensures the dressing remains in place.   |
|   | tions / Precautions  |
| 0 | AEROFILM should not be used on infected wounds unless under the supervision of a healthcare practitioner.  |
| 0 | As with all adhesive products, apply and remove carefully from sensitive or fragile skin.  |



# AEROFILM WATERPROOF FILM DRESSING 6 X 7CM

AFW0063



This product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Waterproof Film Dressing