



ADDITIONAL INFORMATION

Weight	0.03 kg
Dimensions	19 × 9.5 × 15 cm
Carton Qty	96
Pricing Unit	Box/3
Size	6cm x 7cm
Min Order Qty	1
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC Code	42311527
Tax Schedule	Non-Taxable
Item Length	7
Item Width	6
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36



PRODUCT DESCRIPTION

AeroFilm is a sterile, waterproof, ultra-thin transparent film dressing engineered for dependable protection on clean skin, superficial wounds and I.V. devices. The breathable film forms an effective barrier to bacteria and viruses while permitting moisture vapour exchange to help maintain a cleaner environment around the site. Its flexible construction conforms smoothly to contours and sits discreetly under clothing, supporting comfortable wear throughout daily activity in clinical, workplace and at-home settings.



AEROFILM WATERPROOF FILM DRESSING 6 X 7CM BOX/3

AFW0063

SupaBond hypoallergenic adhesive is designed for strong, consistent hold that resists edge lift from movement and moisture. The latex-free formulation supports sensitive skin and removes cleanly when it is time to change. Because the film is crystal clear, the application site can be observed continuously without disturbing the area, helping reduce unnecessary changeovers and supporting good practice monitoring for I.V. insertion sites and minor wound coverage.

An integrated frame applicator simplifies handling and placement for crease-free results with minimal contact to the adhesive surface, even on curved or awkward body contours. The low-profile, ultra-thin film flexes with the body to reduce catching on clothing and helps provide a smooth, low-friction outer surface. Each dressing is sterile and single-use to support hygiene requirements. AeroFilm provides reliable coverage for minor burns, blisters and grazes, and is equally suited to protecting post-procedure areas where a waterproof, breathable barrier is desired.

- **Waterproof barrier:** protects against water and body fluids while remaining breathable for comfort.
- **Secure, skin-friendly hold:** SupaBond hypoallergenic adhesive resists edge lift and is latex-free.
- **Clear visibility:** transparent film enables site checks without removal.
- **Easy placement:** integrated frame applicator aids smooth, crease-free application.
- **Conformable comfort:** ultra-thin film moves with the body and wears discreetly under clothing.

The **6 cm x 7 cm size** offers focused coverage with a low-profile fit on wrists, hands and small dressings. It conforms smoothly to curved areas for discreet protection while allowing quick visual checks.

Packaged in a **box of 3** to support compact first aid kits, travel packs and low-volume clinical use, helping reduce waste while keeping shower-friendly protection on hand.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Waterproof Film Dressing