



ADDITIONAL INFORMATION

Weight	0.0348 kg
Dimensions	2.5 × 9.5 × 13 cm
Carton Qty	120
Pricing Unit	Roll/1
Size	10cm x 1M
Min Order Qty	1
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC Code	42311527
Tax Schedule	Non-Taxable
Item Length	100
Item Width	10
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36



PRODUCT DESCRIPTION

AeroFilm is a sterile, waterproof, ultra-thin transparent film roll designed to protect I.V. sites, clean skin and dressed wounds. The roll format lets teams cut the exact length and width required for joints and awkward areas, reducing waste while ensuring a smooth, custom fit. The film forms a breathable bacterial and viral barrier to outside contaminants and is hospital grade, hypoallergenic and latex free.

SupaBond adhesive and an easy frame-assisted application provide a crease-free seal that conforms to movement. Once applied, the transparent film makes site checks simple without disturbing the area, supporting timely assessments and reducing unnecessary changeovers. The durable, skin-friendly adhesive maintains reliable fixation during daily activity and showering.

Designed for clinical settings and well suited to workplace first aid and home care, the roll secures primary dressings, protects minor burns and blisters, and shields healing skin from friction and moisture. It is sterile and packaged for single use, adding confidence wherever dependable coverage is required.

- **Waterproof, sterile barrier:** protects against water, bacteria and viruses.
- **Cut-to-size roll:** conforms to elbows, knees and other awkward areas.
- **Secure adhesion:** frame-assisted application helps resist edge lift.
- **Transparent and breathable:** enables easy site monitoring and comfort.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.