



## ADDITIONAL INFORMATION

<b>Weight</b>	0.21 kg
<b>Dimensions</b>	25.2 × 5 × 13 cm
<b>Carton Qty</b>	20
<b>Pricing Unit</b>	Box/20
<b>Size</b>	10cm x 20cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC</b>	42311546
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	20
<b>Item Width</b>	10
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	36

## PRODUCT DESCRIPTION

AEROFIX dressings combine a SupaSorb pad with a soft and conformable fixative layer for the simple and effective management of larger wounds, abrasions and sutured wounds.

Non Woven Backing  
comfortable and conforming  
soft non woven backing

Supabond™  
Hypo-allergenic adhesive  
exclusively engineered,  
for low irritation

wound  
exudate pad  
super absorption,  
reducing skin maceration

Non-stick barrier  
protects wound with a non-adherent barrier,  
minimising trauma on removal

## **SupaSorb Wound Pad**

- An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

## **SupaBond adhesive**

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

## **Soft and Conformable**

- A non-woven backing is smooth and soft, and adheres to body contours using SupaBond adhesion.

## **Patient Comfort**

- The SupaSorb non-adherent wound pad , combined with a soft and conforming non-woven backing takes patient comfort to a new level, minimising trauma and assisting recovery.

**Contraindications / Precautions**

- AEROFIX should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

**This product is not a substitute for direct medical attention from a healthcare practitioner.**

Aerofix non-woven Island Dressing