



ADDITIONAL INFORMATION

Weight	0.21 kg
Dimensions	25.2 × 5 × 13 cm
Carton Qty	20
Pricing Unit	Box/20
Size	10cm x 20cm
Min Order Qty	1
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC Code	42311546
Tax Schedule	Taxable
Item Length	20
Item Width	10
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36
Warranty Period (Months)	12



PRODUCT DESCRIPTION

AEROFIX Non-Woven Island Dressing is a sterile, latex-free adhesive dressing designed for simple, effective coverage of larger wounds,

abrasions and sutured sites. A soft, conformable non-woven backing helps the dressing sit smoothly over body contours while maintaining reliable fixation. The island design focuses coverage where it is needed and keeps edges low-profile for comfortable wear under clothing and during movement.

At the centre is the SupaSorb absorbent pad engineered to draw blood and exudate into internal pockets. This high-capacity pad supports a cleaner wound environment and reduces unnecessary changeovers. A non-adherent contact layer sits against the skin to minimise sticking, so removal is gentler and less disruptive to fragile or healing tissue.

Fixation is delivered by SupaBond adhesive, developed for strong, uniform adhesion that helps prevent lift and edge curl. The breathable, water-repellent non-woven barrier enhances patient comfort while providing dependable protection in everyday and clinical settings. Individually sterile and intended for single use, AEROFIX offers clinicians, workplaces and home users a dependable choice for routine wound care.

- **Strong hold:** SupaBond adhesive helps prevent edge lift and displacement.
- **High absorbency:** SupaSorb pad manages exudate to keep the site cleaner.
- **Gentle removal:** non-adherent layer reduces trauma at dressing changes.
- **Conformable comfort:** soft non-woven backing contours to the body.
- **EO-sterile:** suitable for broad clinical and first aid use.

The **9 x 20 cm size** delivers extended coverage for larger post-procedure sites and abrasions while remaining comfortable on curved areas, supporting secure fixation and easier dressing management.

Supplied as a **single sterile dressing**, this option suits topping up kits or trial use where only one island dressing is needed. The individual unit format supports targeted ordering and tailored patient care.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.