



## ADDITIONAL INFORMATION

<b>Weight</b>	0.18 kg
<b>Dimensions</b>	14.7 × 12 × 10 cm
<b>Carton Qty</b>	10
<b>Pricing Unit</b>	Box/50
<b>Size</b>	6cm x 8cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311546
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	8
<b>Item Width</b>	6
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	36
<b>Warranty Period (Months)</b>	12



## PRODUCT DESCRIPTION

AEROFIX Non-Woven Island Dressing is a sterile, latex-free adhesive dressing designed for simple, effective coverage of larger wounds,

abrasions and sutured sites. A soft, conformable non-woven backing helps the dressing sit smoothly over body contours while maintaining reliable fixation. The island design focuses coverage where it is needed and keeps edges low-profile for comfortable wear under clothing and during movement.

At the centre is the SupaSorb absorbent pad engineered to draw blood and exudate into internal pockets. This high-capacity pad supports a cleaner wound environment and reduces unnecessary changeovers. A non-adherent contact layer sits against the skin to minimise sticking, so removal is gentler and less disruptive to fragile or healing tissue.

Fixation is delivered by SupaBond adhesive, developed for strong, uniform adhesion that helps prevent lift and edge curl. The breathable, water-repellent non-woven barrier enhances patient comfort while providing dependable protection in everyday and clinical settings. Individually sterile and intended for single use, AEROFIX offers clinicians, workplaces and home users a dependable choice for routine wound care.

- **Strong hold:** SupaBond adhesive helps prevent edge lift and displacement.
- **High absorbency:** SupaSorb pad manages exudate to keep the site cleaner.
- **Gentle removal:** non-adherent layer reduces trauma at dressing changes.
- **Conformable comfort:** soft non-woven backing contours to the body.
- **EO-sterile:** suitable for broad clinical and first aid use.

The **Size: 6 × 8 cm** option centres a generous absorbent pad within a compact footprint, ideal for small post-procedure sites and everyday cuts where comfortable coverage and secure fixation are required.

Supplied in a **box of 50** EO sterile dressings, this option supports efficient clinical supply and workplace kit stocking by reducing the frequency of reordering and mid shift restocking, while maintaining consistent availability for routine dressing changes.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.