

AEROFIX NON-WOVEN ISLAND DRESSING 9 X 10CM BOX/50

AFID8S





ADDITIONAL INFORMATION

Weight	0.3 kg
Dimensions	15.2 × 12 × 13 cm
Carton Qty	10
Pricing Unit	Box/50
Size	8cm x 10cm
Min Order Qty	1
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	10
Item Width	8
Item Height	12
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36

PRODUCT DESCRIPTION

AEROFIX dressings combine a SupaSorb pad with a soft and conformable fixative layer for the simple and effective management of larger wounds, abrasions and sutured wounds.



Aeroplast™ Non Woven Backing comfortable and conforming soft non woven backing

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation

wound exudate pad super absorption, reducing skin maceration

Non-stick barrier protects wound with a non-adherent barrier, minimising trauma on removal

SupaSorb Wound Pad

 An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

SupaBond adhesive

• SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

Soft and Conformable

• A non-woven backing is smooth and soft, and adheres to body contours using SupaBond adhesion.

Patient Comfort

• The SupaSorb non-adherent wound pad, combined with a soft and conforming non-woven backing takes patient comfort to a new level, minimising trauma and assisting recovery.

Contraindications / Precautions

- AEROFIX should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

AFID8S



AEROFIX NON-WOVEN ISLAND DRESSING 9 X 10CM BOX/50

AFID8S

AeroFix Non-woven Island Dressings 9cm x 10cm Box/50