



ADDITIONAL INFORMATION

Weight	0.023 kg
Dimensions	3.5 × 3.5 × 10 cm
Carton Qty	108
Pricing Unit	Unit/1
Size	10cm x 1M
Min Order Qty	6
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
Tax Schedule	Taxable
Item Length	100
Item Width	10
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60



PRODUCT DESCRIPTION

The AeroFix Non-Woven Fixation Tape represents a pinnacle in medical adhesive technology, offering unparalleled comfort, flexibility, and secure adhesion. This product is expertly crafted to meet the diverse needs of healthcare professionals and patients, ensuring optimal performance across various medical environments. The tape's non-woven fabric is designed for exceptional conformability, making it an ideal choice for securing dressings on joints and other challenging body contours. Its hypoallergenic adhesive is gentle on the skin, making it suitable for sensitive areas, including facial applications.



AEROFIX NON-WOVEN FIXATION TAPE 10CM X 1M

AFR101

Featuring the innovative SupaBond Adhesive Technology, this fixation tape guarantees a strong and reliable hold. Dressings remain securely in place without the risk of lift-up or edge curling, providing peace of mind to both healthcare providers and patients. The AeroFix tape is also designed to promote skin health, thanks to its optimal breathability. The fabric allows air and moisture to pass through, reducing the risk of skin irritation and infections.

Ease of use is a hallmark of the AeroFix Non-Woven Fixation Tape. The application process is straightforward: wash hands before and after application, apply to clean and dry skin without stretching, and press the edges down firmly. The tape should be changed daily to maintain hygiene and effectiveness. Its versatility makes it suitable for a wide range of applications, from securing dressings on sensitive or fragile skin to acting as an underwrap tape for added protection.

Manufactured to the highest standards, this product is CE and UKCA certified, ensuring compliance with international safety and quality regulations. It is also latex-free, making it a safe choice for individuals with latex allergies. The AeroFix Non-Woven Fixation Tape is available in various sizes, including a convenient 10cm x 1M roll, providing ample coverage for different medical needs. Whether used in hospitals, clinics, or at home, this tape offers a reliable solution for dressing retention, enhancing the overall patient care experience.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Underwrap Tape 10cm x 1M