



ADDITIONAL INFORMATION

Weight	0.013 kg
Dimensions	3.5 × 3.5 × 5.5 cm
Carton Qty	168
Pricing Unit	Unit/1
Size	5cm x 1M
Min Order Qty	6
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
Tax Schedule	Taxable
Item Length	100
Item Width	5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

Our specially designed fixation tape offers comfort and flexibility for every contour of the body, making it the ideal choice for dressings, sensitive areas, and as an underwrap tape.

SupaBond Adhesive Technology: Our unique SupaBond adhesive guarantees exceptional adhesion, ensuring secure placement without any lift-up or edge curling.

High Conformability: AeroFix UltraFlex smoothly adheres to body contours, ensuring optimal fixation even over joints and tricky body

areas.

Optimal Breathability: The non-woven fabric promotes skin health by allowing it to breathe, thereby reducing the risk of infections.

Gentle On The Skin: Designed to be hypo-allergenic and soft, it is perfect for sensitive, fragile, or facial skin areas. However, apply and remove carefully from such skin to ensure utmost safety.

Contraindications / Precautions

- AeroFix should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Underwrap Tape 5cm x 1M