



**ADDITIONAL INFORMATION**

<b>Weight</b>	0.013 kg
<b>Dimensions</b>	3.5 × 3.5 × 5.5 cm
<b>Carton Qty</b>	168
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	5cm x 1M
<b>Min Order Qty</b>	6
<b>ARTG ID</b>	236937
<b>GMDN Code</b>	34830
<b>Medical Device Class</b>	Class 1
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	100
<b>Item Width</b>	5
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Barcode (GTIN/EAN/UPC)</b>	9341394028777

**PRODUCT DESCRIPTION**

Our specially designed fixation tape offers comfort and flexibility for every contour of the body, making it the ideal choice for dressings, sensitive areas, and as an underwrap tape.

**SupaBond Adhesive Technology:** Our unique SupaBond adhesive guarantees exceptional adhesion, ensuring secure placement without any lift-up or edge curling.

**High Conformability:** AeroFix UltraFlex smoothly adheres to body contours, ensuring optimal fixation even over joints and tricky body areas.

**Optimal Breathability:** The non-woven fabric promotes skin health by allowing it to breathe, thereby reducing the risk of infections.

**Gentle On The Skin:** Designed to be hypo-allergenic and soft, it is perfect for sensitive, fragile, or facial skin areas. However, apply and remove carefully from such skin to ensure utmost safety.

## **Contraindications / Precautions**

- AeroFix should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

**This product is not a substitute for direct medical attention from a healthcare practitioner.**

Underwrap Tape 5cm x 1M