



ADDITIONAL INFORMATION

Weight	0.02 kg
Dimensions	5.2 × 5.1 × 6.5 cm
Carton Qty	120
Pricing Unit	Box/1
Size	6cm x 5M
Min Order Qty	1
ARTG ID	236937
GMDN Code	34830
Medical Device Class	Class 1
Tax Schedule	Taxable
Item Length	
Item Width	
Item Measurement Unit	m
Location	AU-Cameron Dr
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AeroForm Cohesive Bandage is designed to provide secure dressing retention with a gentle, non-slip hold. Made from lightweight gauze, it offers a balance of flexibility and support, making it suitable for use on limbs and contoured areas. The bandage forms a stable layer that helps keep primary dressings in place during everyday movement, supporting consistent coverage without relying on clips, tapes or pins.

The cohesive structure allows the bandage to adhere to itself rather than to skin or hair. This helps reduce discomfort during

adjustment or removal and makes the product a practical option for people with sensitive or fragile skin. The construction supports repeated wrapping to build up layers where additional compression or stability is required, such as around joints or over padded dressings.

Breathable gauze fabric helps air circulate around the covered area, supporting comfort during extended wear. The material is latex free, reducing the risk of reactions for those with latex sensitivity. The flexible nature of the bandage allows it to move with the body, helping to maintain secure retention while supporting a natural range of motion in knees, elbows, wrists and ankles.

Across the range, AeroForm Cohesive Bandage is suitable for dressing retention in workplaces, clinics, sports settings and home first aid. It is appropriate for situations where a firm yet gentle hold is required, such as securing absorbent pads, gauze swabs or protective coverings over minor injuries. The cohesive design also makes it useful when frequent inspection or redressing is expected, as the bandage can be reapplied with minimal disruption.

- **Cohesive hold:** adheres to itself, not to skin or hair, for more comfortable removal.
- **Gentle on sensitive skin:** latex-free gauze construction helps reduce irritation.
- **Breathable fabric:** lightweight material promotes airflow for extended wear comfort.
- **Flexible support:** conforms around limbs and joints while retaining dressings securely.
- **Layerable design:** can be wrapped in multiple layers to tailor support and retention.

The **Size – 6cm x 2m** format offers balanced coverage for limbs and medium areas, allowing comfortable wrapping around joints while maintaining an easily managed width for precise dressing retention.

Supplied in a **single roll** this option is ideal for adding a cohesive bandage to a first aid kit or keeping one readily available for everyday use. A single roll provides the same breathable, latex-free gauze and cohesive self-grip performance, suitable for individual dressing changes or treating one person as needed.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.