



ADDITIONAL INFORMATION

Weight	0.06 kg
Dimensions	21 × 20 × 3.5 cm
Carton Qty	30
Pricing Unit	Wrap/12
Size	10cm x 4M
Min Order Qty	1
ARTG ID	236936
GMDN Code	48126
Medical Device Class	Class 1
UNSPSC	4231 1500
Tax Schedule	Taxable
Item Length	400
Item Width	10
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60

PRODUCT DESCRIPTION

The AeroForm Conforming bandage is a quality cost-effective bandage with excellent conforming attributes. Our Heavyweight and Premium range offers greater support and compression and are recommended for sprains, strains and snake bites.

Conforming Specifications

AeroForm Conforming Bandages are of a standard conforming nature which is durable and suitable for dressing retention.

Easy to Open Wrapper



AEROFORM CONFORMING BANDAGE 10CM X 4M WRAP/12

AF100

AeroForm Conforming Bandages have a easy open tab built into the wrapper of the bandage.

Conforming Bandage 10cm