



## ADDITIONAL INFORMATION

|                                 |                  |
|---------------------------------|------------------|
| <b>Weight</b>                   | 0.020833 kg      |
| <b>Dimensions</b>               | 30 × 21 × 3.5 cm |
| <b>Carton Qty</b>               | 20               |
| <b>Pricing Unit</b>             | Wrap/12          |
| <b>Size</b>                     | 15cm x 4M        |
| <b>Min Order Qty</b>            | 1                |
| <b>ARTG ID</b>                  | 236936           |
| <b>GMDN Code</b>                | 48126            |
| <b>Medical Device Class</b>     | Class 1          |
| <b>UNSPSC Code</b>              | 42311500         |
| <b>Tax Schedule</b>             | Taxable          |
| <b>Item Length</b>              | 400              |
| <b>Item Width</b>               | 15               |
| <b>Item Measurement Unit</b>    | cm               |
| <b>Is Medical Device</b>        | Yes              |
| <b>Is Serialized Item</b>       | No               |
| <b>Is Dangerous Good</b>        | No               |
| <b>Shelf Life (Months)</b>      | 60               |
| <b>Warranty Period (Months)</b> | 0                |



## PRODUCT DESCRIPTION

AeroForm Conforming Bandage is a lightweight, elastic gauze bandage designed for reliable dressing retention and gentle, even

compression without constriction. The soft, breathable weave conforms closely to body contours so pads and dressings stay securely in place on wrists, ankles and other curved or moving areas. A clean, non-sterile, single-use format supports everyday first aid, workplace incidents and clinical dressing changes where a conforming hold is needed.

The bandage stretches smoothly to create a low-profile wrap that is easy to apply and comfortable to wear. The elastic construction helps maintain consistent pressure while accommodating natural swelling and movement. Latex-free composition and skin-friendly edges support use on a wide range of patients, including sensitive users. Packaging features an easy-open tab for swift access when time matters.

Engineered for dependable performance across care settings, this conforming bandage is suited to retaining primary wound dressings, securing splints and padding, and providing light compression during recovery. It is supplied in flow-wrapped singles with shrink-wrapped inners for efficient inventory control. Ideal for stocking in first aid kits, treatment rooms, ambulances and event medical stations.

- **Conforming elasticity:** wraps smoothly to body contours for a secure, comfortable hold
- **Lightweight and breathable:** soft gauze construction supports comfort during wear
- **Latex-free and hypoallergenic:** suitable for sensitive skin policies
- **Gentle compression:** provides support without constriction for everyday use
- **Easy-open packaging:** fast access with flow wrap easy-open tab
- **Single-use, non-sterile:** appropriate for routine dressing retention
- **Versatile application:** suitable for first aid, workplace and clinical settings

The **Size 15 cm x 4 m** option offers broad coverage that is ideal for wrapping larger areas such as lower legs and thighs, delivering gentle, even support while helping to keep primary dressings securely in place.

This size offers a **15 cm width** and **4 m stretched length**, providing generous coverage for limbs and joints. The broader width supports secure dressing retention on larger areas while maintaining comfortable, light compression.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Conforming Bandage 15cm