



ADDITIONAL INFORMATION

Weight	0.008333 kg
Dimensions	21 × 15 × 3.5 cm
Carton Qty	40
Pricing Unit	Wrap/12
Size	7.5cm x 4M
Min Order Qty	1
ARTG ID	236936
GMDN Code	48126
Medical Device Class	Class 1
UNSPSC Code	42311500
Tax Schedule	Taxable
Item Length	400
Item Width	7.5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60
Warranty Period (Months)	0



PRODUCT DESCRIPTION

AeroForm Conforming Bandage is a lightweight, elastic gauze bandage designed for reliable dressing retention and gentle, even

compression without constriction. The soft, breathable weave conforms closely to body contours so pads and dressings stay securely in place on wrists, ankles and other curved or moving areas. A clean, non-sterile, single-use format supports everyday first aid, workplace incidents and clinical dressing changes where a conforming hold is needed.

The bandage stretches smoothly to create a low-profile wrap that is easy to apply and comfortable to wear. The elastic construction helps maintain consistent pressure while accommodating natural swelling and movement. Latex-free composition and skin-friendly edges support use on a wide range of patients, including sensitive users. Packaging features an easy-open tab for swift access when time matters.

Engineered for dependable performance across care settings, this conforming bandage is suited to retaining primary wound dressings, securing splints and padding, and providing light compression during recovery. It is supplied in flow-wrapped singles with shrink-wrapped inners for efficient inventory control. Ideal for stocking in first aid kits, treatment rooms, ambulances and event medical stations.

- **Conforming elasticity:** wraps smoothly to body contours for a secure, comfortable hold
- **Lightweight and breathable:** soft gauze construction supports comfort during wear
- **Latex-free and hypoallergenic:** suitable for sensitive skin policies
- **Gentle compression:** provides support without constriction for everyday use
- **Easy-open packaging:** fast access with flow wrap easy-open tab
- **Single-use, non-sterile:** appropriate for routine dressing retention
- **Versatile application:** suitable for first aid, workplace and clinical settings

The **7.5 cm x 4 m size** offers broader coverage for medium areas of the body, creating a comfortable wrap that conforms smoothly while providing light, even support.

This bandage is supplied as 12 rolls 7.5 cm x 4 m, giving balanced coverage for medium areas such as wrists, forearms and ankles. The length allows multiple overlapping turns for secure dressing retention and light, even compression while remaining easy to handle in first aid situations.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Conforming Bandage 7.5cm