


### **AEROFORM SHORT SNAKE BITE BANDAGE WITH INDICATOR 10CM X 4.5M WRAP/12**

AFSB10





#### **ADDITIONAL INFORMATION**

Weight 0.52 kg

**Dimensions**  $22 \times 4 \times 26$  cm

**Carton Qty** 20

**Pricing Unit** Wrap/12

Size 10cm x 4.5M

**Min Order Qty** 

**ARTG ID** 236934 **GMDN** Code 10279 **Medical Device Class** Class 1

UNSPSC 42311506

**Tax Schedule** Taxable

**Item Length** 450

**Item Width** 10

Item Measurement Unit cm **Is Medical Device** 

Yes

Is Serialized Item No

Is Dangerous Good No

#### **PRODUCT DESCRIPTION**

The AeroForm Premium Snake Bite Bandage with Continuous Indicators is a Premium Heavyweight Conforming Bandage. Perfect for use with the Pressure Immobilisation Technique, which is used for the treatment of all Australian Snake and Funnel Webbed Spider bites, Blue Ringed Octopus and Cone Shell stings.

This bandage is of sufficient length to bandage a large adult leg.



## AEROFORM SHORT SNAKE BITE BANDAGE WITH INDICATOR 10CM X 4.5M WRAP/12

AFSB10

○ 10cm x 4.5M Stretched
Benefits
Wider and longer than most crepe or elastic bandages.
<ul> <li>Continuous indicators ensure you have the correct pressure over the entire limb.</li> </ul>
<ul> <li>Suitable for the treatment of all Australian Snake and Funnel Webbed Spider bites, Blue Ringed Octopus and Cone Shell stings.</li> </ul>
Indicators show when the correct level of compression is achieved
The unique stitching technology ensures the correct pressure is evenly applied over the entire limb.
This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



# AEROFORM SHORT SNAKE BITE BANDAGE WITH INDICATOR 10CM X 4.5M WRAP/12

AFSB10

Short Snake Bite Bandage 10cm