



ADDITIONAL INFORMATION

Weight	0.66 kg
Dimensions	20 × 14 × 11.5 cm
Carton Qty	12
Pricing Unit	Box/6
Size	10cm x 10.5M
Min Order Qty	1
ARTG ID	515586
GMDN Code	10284
Medical Device Class	Class 1
UNSPSC Code	42311506
Tax Schedule	Taxable
Item Length	1050
Item Width	10
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AeroForm Snake Bite Bandage with Indicator is a premium, heavyweight conforming elastic bandage designed to support pressure immobilisation for venomous bites and stings. Printed tension indicators provide a clear visual cue to help reach effective compression along the limb, offering reassurance when every second counts in the outdoors, on worksites, and at home. The indicator pattern is easy to recognise, helping users confirm when the correct extension has been reached for stable support during immobilisation.

The durable elastic fabric conforms smoothly to body contours while maintaining firm, even support with movement. The robust

weave resists bunching as you apply tension, helping to create a smooth, low-profile wrap with reliable overlap and minimal gaps. Clear markings assist responders and trainees alike to achieve consistent technique and repeatable results under pressure. Non-sterile and single use, the bandage is a dependable inclusion for first aid kits in vehicles, workplaces, community organisations, and household preparedness.

Built for practical performance in real-world conditions, AeroForm is purpose-made for emergency environments where clarity and control matter. It supports uniform pressure across the covered area and pairs well with limb immobilisation methods until professional medical care is available. The construction balances a strong hold with comfortable conformity, delivering robust compression without adding unnecessary bulk—ideal for outdoor adventures, remote travel, industrial settings, and community events.

- **Visual tension guide:** indicator prints show when effective compression is achieved.
- **Heavyweight elastic build:** firm, even support that conforms to the limb.
- **Single use, non-sterile:** practical for well-stocked first aid kits.
- **Smooth, stable wrap:** robust weave resists bunching for consistent overlap.
- **Confidence under pressure:** helps users apply reliable compression quickly.

Size – Long 10 cm x 10.5 m: designed for comprehensive limb coverage with heavyweight elastic that conforms smoothly for stable, even compression.

Supplied in a **box of 6**, this option makes it easy to stock multiple snake bite bandages in first aid kits for workplaces, vehicles, outdoor programs, and community groups, ensuring several wraps are immediately available when needed.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Indicator Snake Bite Bandage Box/6