



## ADDITIONAL INFORMATION

<b>Weight</b>	0.019 kg
<b>Dimensions</b>	11.5 × 4.9 × 0.5 cm
<b>Carton Qty</b>	1008
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	11cm
<b>Min Order Qty</b>	12
<b>ARTG ID</b>	236810
<b>GMDN Code</b>	13481
<b>Medical Device Class</b>	Class 1
<b>UNSPSC Code</b>	42291614
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	11
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



## PRODUCT DESCRIPTION

AEROINSTRUMENTS Stainless Steel Scissors deliver reliable, precise cutting for first aid tasks and clinical support activities. Built from quality stainless steel, the blades align accurately at the pivot to provide a smooth, controlled cut through dressings, gauze, medical tape and light fabrics encountered during first aid. The compact design fits neatly into workplace kits, ambulance bags and home cabinets, making it a practical choice for everyday readiness.

A polished, corrosion-resistant finish supports long service life in demanding environments, while rounded finger and thumb rings offer familiar control for confident handling. The instrument is supplied non-sterile and is intended for general first aid use, including trimming of bandages and clothing as required during assessment and treatment—this product is not intended for surgical use.

Consistent build quality helps maintain dependable performance between replenishments across the range.

Trusted by organisations assembling compliant first aid sets, these scissors balance portability and durability without unnecessary bulk. Assorted sizes and specifications are available across the range to meet different kit layouts and cutting tasks, supporting standardisation for workplaces, education, and community responders.

- **Durable stainless steel:** corrosion-resistant construction for long-lasting use.
- **Controlled cutting:** precisely ground blades for smooth, accurate performance.
- **Comfortable grip:** classic finger and thumb rings for stable handling.
- **First aid ready:** suitable for trimming dressings, tape and light fabric; non-sterile, not for surgical use.

**Plastic Handle** option offers a lightweight, comfortable grip with classic rings for steady control during precise trimming in first aid applications.

The **11 cm length** offers compact control that fits neatly into first aid kits and pockets while providing precise handling for trimming dressings, gauze and tape during quick assessments.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Scissors - S/S - Plastic Handle