



ADDITIONAL INFORMATION

Weight	0.02 kg
Dimensions	9.5 × 0.5 × 4.4 cm
Carton Qty	400
Pricing Unit	Unit/1
Size	9cm
Min Order Qty	50
ARTG ID	236810
GMDN Code	13481
Medical Device Class	Class 1
UNSPSC Code	42291614
Tax Schedule	Taxable
Item Length	9
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AEROINSTRUMENTS Stainless Steel Scissors deliver reliable, precise cutting for first aid tasks and clinical support activities. Built from quality stainless steel, the blades align accurately at the pivot to provide a smooth, controlled cut through dressings, gauze, medical tape and light fabrics encountered during first aid. The compact design fits neatly into workplace kits, ambulance bags and home cabinets, making it a practical choice for everyday readiness.

A polished, corrosion-resistant finish supports long service life in demanding environments, while rounded finger and thumb rings offer familiar control for confident handling. The instrument is supplied non-sterile and is intended for general first aid use, including trimming of bandages and clothing as required during assessment and treatment—this product is not intended for surgical use.

Consistent build quality helps maintain dependable performance between replenishments across the range.

Trusted by organisations assembling compliant first aid sets, these scissors balance portability and durability without unnecessary bulk. Assorted sizes and specifications are available across the range to meet different kit layouts and cutting tasks, supporting standardisation for workplaces, education, and community responders.

- **Durable stainless steel:** corrosion-resistant construction for long-lasting use.
- **Controlled cutting:** precisely ground blades for smooth, accurate performance.
- **Comfortable grip:** classic finger and thumb rings for stable handling.
- **First aid ready:** suitable for trimming dressings, tape and light fabric; non-sterile, not for surgical use.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Scissors - S/S - Sharp/Sharp