

AEROMASK CHILDREN'S SURGICAL MASK BOX/50

ASMC50



ADDITIONAL INFORMATION

Weight 0.0032 kg

Dimensions $10 \times 15.5 \times 10.3$ cm

Carton Qty 56

Pricing Unit Box/50

Size Small – 14.5cm x 9cm

Min Order Qty 1

ARTG ID 339645
GMDN Code 35177
Medical Device Class 1

UNSPSC 42131713

Tax Schedule Taxable

Item Length 14.5

Item Width 9

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

Children's Surgical Mask

The AeroMask ASMC50 face mask is a premium quality 3ply surgical mask with Level 2 protection rating. The masks are high quality with earloops and come with nose clips to ensure secure fitting.

Size: Small - 14.5cm x 9cm

3 Ply

Level 2 rating BFE ≥ 98% Elastic earloops Premium quality Disposable Latex free

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the



AEROMASK CHILDREN'S SURGICAL MASK BOX/50

ASMC50

	٠					_	
М	1	rΔ	cti	r	10	tor	use.
u	ı	ıc	СU	OI.	ıs	101	use.

AEROMASK Children's Surgical Masks