



**ADDITIONAL INFORMATION**

<b>Weight</b>	0.0032 kg
<b>Dimensions</b>	10 × 15.5 × 10.3 cm
<b>Carton Qty</b>	56
<b>Pricing Unit</b>	Box/50
<b>Size</b>	Small - 14.5cm x 9cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	339645
<b>GMDN Code</b>	35177
<b>Medical Device Class</b>	Class 1
<b>UNSPSC</b>	42131713
<b>Tax Schedule</b>	Taxable
<b>Location</b>	AU-Warehouse
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	F
<b>Groups of Qty</b>	1
<b>Item Length</b>	14.5
<b>Item Width</b>	9
<b>Item Measurement Unit</b>	cm

**PRODUCT DESCRIPTION**

Children's Surgical Mask

The AeroMask ASMC50 face mask is a premium quality 3ply surgical mask with Level 2 protection rating. The masks are high quality with earloops and come with nose clips to ensure secure fitting.

**Size: Small - 14.5cm x 9cm**

**3 Ply**

Level 2 rating



# AEROMASK CHILDREN'S SURGICAL MASK BOX/50

ASMC50

BFE  $\geq$  98%

Elastic earloops

Premium quality

Disposable

Latex free

AEROMASK Children's Surgical Masks