



ADDITIONAL INFORMATION

| | |
|-----------------------------|----------------|
| Weight | 0.056 kg |
| Dimensions | 4 × 13 × 10 cm |
| Carton Qty | 100 |
| Pricing Unit | Unit/1 |
| Min Order Qty | 1 |
| ARTG ID | 236803 |
| GMDN Code | 36066 |
| Medical Device Class | Class 1 |
| UNSPSC Code | 42272303 |
| Tax Schedule | Non-Taxable |
| Is Medical Device | Yes |
| Is Serialized Item | No |
| Is Dangerous Good | No |
| Shelf Life (Months) | 36 |



PRODUCT DESCRIPTION

AEROMASK provides a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

CONTENTS

1. Folded cushion mask with Oxygen Port
2. One-way valve

3. Disposable filter
4. Elastic head strap

All materials are Latex Free

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Pocket Mask in Bag