



ADDITIONAL INFORMATION

Weight 0.056 kg

Dimensions $4 \times 13 \times 10 \text{ cm}$

Carton Qty 100

Pricing Unit Unit/1

Min Order Qty 1

ARTG ID 236803 **GMDN Code** 36066

Medical Device Class 1

UNSPSC 42272303

Tax Schedule Non-Taxable

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No
Shelf Life (Months) 36

PRODUCT DESCRIPTION

AEROMASK provides a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

CONTENTS

- 1. Folded cushion mask with Oxygen Port
- 2. One-way valve
- 3. Disposable filter
- 4. Elastic head strap



AEROSHIELD CPR MASK IN CLEAR BAG

AM01B

Αll	materials	are	Latex	Free

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Pocket Mask in Bag