



ADDITIONAL INFORMATION

Weight	0.056 kg
Dimensions	4 × 13 × 10 cm
Carton Qty	100
Pricing Unit	Unit/1
Min Order Qty	1
ARTG ID	236803
GMDN Code	36066
Medical Device Class	Class 1
UNSPSC	42272303
Tax Schedule	Non-Taxable
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394008045

PRODUCT DESCRIPTION

AEROMASK provides a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

CONTENTS

1. Folded cushion mask with Oxygen Port
2. One-way valve
3. Disposable filter
4. Elastic head strap

All materials are Latex Free



AEROSHIELD CPR MASK IN CLEAR BAG

AM01B

Pocket Mask in Bag