



ADDITIONAL INFORMATION

Weight	0.043 kg
Dimensions	10 × 20 × 2.5 cm
Carton Qty	120
Pricing Unit	Box/10
Size	One Size
Min Order Qty	1
ARTG ID	339645
GMDN Code	35177
Medical Device Class	Class 1
UNSPSC	42131713
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

ARTG No. 339645

AEROMASK provides a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

Surgical Mask

The AeroMask SM01 face mask is a premium quality 3ply surgical mask with Level 2 protection rating. The masks are high quality with earloops and come with nose clips to ensure secure fitting. The SM01 mask in boxes of 10 is known as the ASM10.

3 Ply

Level 2 rating

BFE ≥ 98%

Elastic earloops

Premium quality



AEROMASK LEVEL 2 SURGICAL MASK BOX/10

ASM10

Disposable

Latex free

AEROMASK Surgical Masks