



**ADDITIONAL INFORMATION**

<b>Weight</b>	0.043 kg
<b>Dimensions</b>	10 × 20 × 2.5 cm
<b>Carton Qty</b>	120
<b>Pricing Unit</b>	Box/10
<b>Size</b>	One Size
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	339645
<b>GMDN Code</b>	35177
<b>Medical Device Class</b>	Class 1
<b>UNSPSC Code</b>	42131713
<b>Tax Schedule</b>	Taxable
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



**PRODUCT DESCRIPTION**

ARTG No. 339645

AEROMASK provides a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

**Surgical Mask**

The AeroMask SM01 face mask is a premium quality 3ply surgical mask with Level 2 protection rating. The masks are high quality with earloops and come with nose clips to ensure secure fitting. The SM01 mask in boxes of 10 is known as the ASM10.



# AEROMASK LEVEL 2 SURGICAL MASK BOX/10

ASM10

## 3 Ply

Level 2 rating

**BFE  $\geq$  98%**

Elastic earloops

Premium quality

Disposable

Latex free

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

AEROMASK Surgical Masks