



ADDITIONAL INFORMATION

| | |
|-----------------------------|------------------|
| Weight | 0.043 kg |
| Dimensions | 10 × 20 × 2.5 cm |
| Carton Qty | 120 |
| Pricing Unit | Box/10 |
| Size | One Size |
| Min Order Qty | 1 |
| ARTG ID | 339645 |
| GMDN Code | 35177 |
| Medical Device Class | Class 1 |
| UNSPSC | 42131713 |
| Tax Schedule | Taxable |
| Location | AU-Warehouse |
| Is Medical Device | Yes |
| Is Serialized Item | No |

PRODUCT DESCRIPTION

ARTG No. 339645

AEROMASK provides a selection of protective and airway management products for pre-hospital treatment on injured patients, and to assist in the safe application of essential CPR.

Surgical Mask

The AeroMask SM01 face mask is a premium quality 3ply surgical mask with Level 2 protection rating. The masks are high quality with earloops and come with nose clips to ensure secure fitting. The SM01 mask in boxes of 10 is known as the ASM10.

3 Ply

Level 2 rating

BFE ≥ 98%

Elastic earloops



AEROMASK LEVEL 2 SURGICAL MASK BOX/10

ASM10

Premium quality

Disposable

Latex free

AEROMASK Surgical Masks