



## ADDITIONAL INFORMATION

<b>Weight</b>	0.1 kg
<b>Dimensions</b>	8.3 × 12.5 × 10.2 cm
<b>Carton Qty</b>	30
<b>Pricing Unit</b>	Box/50
<b>Size</b>	5cm x 5cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	159437
<b>GMDN Code</b>	34654
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311545
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	5
<b>Item Width</b>	5
<b>Item Height</b>	8
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60
<b>Warranty Period (Months)</b>	0



## PRODUCT DESCRIPTION

AeroPad™ Non-Adherent Dressing is a soft, non-stick wound pad designed to protect delicate skin while offering reliable absorbency. A dual-sided non-adherent contact layer helps minimise trauma during dressing changes, while the cushioned inner core draws fluid away from the site to support comfort. The dressing is sterile and latex free, making it suitable for use across workplaces, schools, clinics and home first aid kits.

The pad sits comfortably beneath securement such as film dressings, paper tapes or conforming bandages, allowing teams to build flexible wound-cover solutions for cuts, grazes and minor burns. Its low-lint, non-woven construction conforms to body contours and is easy to position for effective coverage on hands, joints and other high-movement areas.

Engineered for dependable performance, AeroPad™ helps standardise wound care across busy environments. The combination of a gentle non-adherent surface and an efficient absorbent layer supports practical changeovers and user comfort, making it a smart choice for day-to-day first aid readiness and clinical triage.

- **Non-adherent contact layer:** helps reduce pain and trauma on removal.
- **Absorbent inner core:** manages exudate for drier, more comfortable wear.
- **Sterile and latex free:** suitable for sensitive users and diverse settings.
- **Soft, conformable pad:** low lint design sits neatly under tape or film.

**Size: 5 × 5 cm** offers compact, square coverage for small cuts and grazes while remaining comfortable on hands and joints. The pad layers neatly under tape or film for streamlined first aid.

Supplied in a **box of 50** to simplify restocking for clinics, workplaces and kits, reducing mid-week top-ups and ensuring a consistent supply for everyday minor wound care.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.