



ADDITIONAL INFORMATION

Weight	0.17 kg
Dimensions	10.8 × 12.5 × 12.8 cm
Carton Qty	10
Pricing Unit	Box/50
Size	7.5cm x 5cm
Min Order Qty	1
ARTG ID	159437
GMDN Code	34654
Medical Device Class	Class 1s
UNSPSC Code	42311545
Tax Schedule	Taxable
Item Length	7.5
Item Width	5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60
Warranty Period (Months)	0



PRODUCT DESCRIPTION

AeroPad™ Non-Adherent Dressing is a soft, non-stick wound pad designed to protect delicate skin while offering reliable absorbency. A

dual-sided non-adherent contact layer helps minimise trauma during dressing changes, while the cushioned inner core draws fluid away from the site to support comfort. The dressing is sterile and latex free, making it suitable for use across workplaces, schools, clinics and home first aid kits.

The pad sits comfortably beneath securement such as film dressings, paper tapes or conforming bandages, allowing teams to build flexible wound-cover solutions for cuts, grazes and minor burns. Its low-lint, non-woven construction conforms to body contours and is easy to position for effective coverage on hands, joints and other high-movement areas.

Engineered for dependable performance, AeroPad™ helps standardise wound care across busy environments. The combination of a gentle non-adherent surface and an efficient absorbent layer supports practical changeovers and user comfort, making it a smart choice for day-to-day first aid readiness and clinical triage.

- **Non-adherent contact layer:** helps reduce pain and trauma on removal.
- **Absorbent inner core:** manages exudate for drier, more comfortable wear.
- **Sterile and latex free:** suitable for sensitive users and diverse settings.
- **Soft, conformable pad:** low lint design sits neatly under tape or film.

Size: 5 × 7.5 cm provides versatile coverage for small to medium wounds while remaining comfortable on hands and joints. The low-profile pad fits neatly under film or tape for dependable everyday protection.

Supplied as a **single unit**, ideal for topping up first aid kits or replacing a single dressing. The individual sterile pad reduces waste and keeps replenishment simple for on-the-spot care.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.