



ADDITIONAL INFORMATION

Weight	0.11111 kg
Dimensions	4.1 × 12 × 10 cm
Carton Qty	36
Pricing Unit	Box/100
Size	Assorted
Min Order Qty	12
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60

PRODUCT DESCRIPTION

The AEROPLAST Plastic dressing is a wash proof sheer bandage recommended for all general applications.

SupaSorb™ Wound Pad

- An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

SupaBond™ Adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

Conforming and Air Permeable

- The AeroPlast Plastic Dressing is made of unique materials that bend, stretch and conform to different shapes and forms, as well as being air permeable which allows the skin and wound to breathe as it heals.

Contraindications / Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.



AEROPLAST PLASTIC X-WIDE STRIP 7.5 X 2.5CM BOX/100

AP801