



ADDITIONAL INFORMATION

Weight	0.3 kg
Dimensions	19.5 × 17.5 × 5.5 cm
Carton Qty	16
Pricing Unit	Bag/500
Size	7.2 x 1.9cm
Min Order Qty	4
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	7.2
Item Width	1.9
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394006669

PRODUCT DESCRIPTION

The AEROPLAST Plastic dressing is a wash proof sheer bandage recommended for all general applications.

SupaSorb™ Wound Pad

- o An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the

wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

SupaBond™ Adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

Conforming and Air Permeable

- The AeroPlast Plastic Dressing is made of unique materials that bend, stretch and conform to different shapes and forms, as well as being air permeable which allows the skin and wound to breathe as it heals.

Contraindications / Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.

CONTENTS LISTS

Code	Description	Qty
AP101	AEROPLAST Plastic Standard Strip 7.2 x 1.9cm Box/100	5
RES205255	Clear Resealable Bag 20.5 x 25.5cm (10" x 8")	1

PLASTIC STANDARD STRIP BAG/500



AEROPLAST PLASTIC STANDARD STRIP 7.2 X 1.9CM BAG/500

AP1500