

AEROPLAST PREMIUM DETECTABLE FINGER EXTENSION PLASTERS 16 X 2CM BOX/50

ADDITIONAL INFORMATION

| Weight | 0.1 kg |
|-----------------------|--------------------|
| Dimensions | 16.5 × 4 × 10.5 cm |
| Carton Qty | 72 |
| Pricing Unit | Box/50 |
| Size | 16.0cm x 2.0cm |
| Min Order Qty | 12 |
| ARTG ID | 151434 |
| GMDN Code | 44990 |
| Medical Device Class | Class 1s |
| UNSPSC | 42311546 |
| Tax Schedule | Taxable |
| Item Length | 16 |
| Item Width | 2 |
| Item Measurement Unit | cm |
| Is Medical Device | Yes |
| Is Serialized Item | No |
| Is Dangerous Good | No |
| | |





PRODUCT DESCRIPTION

The Steroplast Sterochef Blue Detectable plaster is a hypo-allergenic, visual and metal detectable plaster specifically designed for the food manufacturing and processing sector.

Visually and Metal Detectable

The Steroplast Sterocheck blue detectable finger extension plaster has multiple detectable attributes. The first being it 's vibrant blue colour, which is the only colour that does not occur naturally in food, making it easily visually detectable. The second feature is the ability to be detected where metal detection equipment is operational. The ultimate in high security.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AEROPLAST PREMIUM DETECTABLE FINGER EXTENSION PLASTERS 16 X 2CM BOX/50

5914

Finger Extension Plaster