



ADDITIONAL INFORMATION

Weight	0.1 kg
Dimensions	16.5 × 4 × 10.5 cm
Carton Qty	72
Pricing Unit	Box/50
Size	16.0cm x 2.0cm
Min Order Qty	12
ARTG ID	151434
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	16
Item Width	2
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The Steroplast SteroChef Blue Detectable plaster is a hypo-allergenic, visual and metal detectable plaster specifically designed for the food manufacturing and processing sector.

Visually and Metal Detectable

The Steroplast Sterocheck blue detectable finger extension plaster has multiple detectable attributes. The first being it 's vibrant blue colour, which is the only colour that does not occur naturally in food, making it easily visually detectable. The second feature is the ability to be detected where metal detection equipment is operational. The ultimate in high security.

Finger Extension Plaster



AEROPLAST PREMIUM DETECTABLE FINGER EXTENSION PLASTERS 16 X 2CM BOX/50

5914