



ADDITIONAL INFORMATION

Weight	0.026 kg
Dimensions	0.5 × 9 × 11 cm
Carton Qty	375
Pricing Unit	Env/20
Size	Assorted
Min Order Qty	15
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394007147

PRODUCT DESCRIPTION

The AeroPlast™ Premium Weight Fabric Adhesive Bandage is a top-tier, durable dressing designed for extended wear and maximum protection. Featuring our SupaBond™ extra-strong adhesive, this bandage is the ideal choice when long-lasting, reliable coverage is needed.

Aeroplast™
Premium Weight Fabric
maximum long wearing durable,
fabric for ultimate wound protection

Supabond™
Hypo-allergenic adhesive
exclusively engineered,
for low irritation and high adhesion

wound
exudate pad
super absorption,
reducing skin maceration

Non-stick barrier
protects wound with a non-adherent barrier,
minimising trauma on removal

Unmatched SupaBond™ Adhesive

Engineered exclusively for AeroPlast™, the SupaBond™ adhesive offers unrivaled adhesion that secures firmly across the application site, preventing lift-up or edge curling. SupaBond™ is latex-free and leaves no sticky residue, addressing common issues found with heavy-duty fabric plasters.

Superior SupaSorb™ Wound Pad The SupaSorb™ wound pad is designed for optimal wound management, effectively absorbing blood and exudates into its internal pockets. This super-absorbent pad reduces the frequency of dressing changes, protecting the wound and keeping it clean and comfortable. The non-adherent film layer minimizes pain and trauma during dressing changes.

Heavy-Duty Fabric Backing

Constructed with conforming and air-permeable heavy-duty fabric, the AEROPLAST™ Premium Fabric dressing is perfect for demanding applications, such as in the manufacturing industry. Its superior adhesion and flexibility ensure the bandage stays in place and conforms to body contours, maintaining wound cleanliness and comfort under the toughest conditions.

The assorted dressings provides a range of popular sized bandages suitable for various injury locations and laceration types. Offering a durable protective layer to create a moist healing environment.

Available as an envelope of 20

Contraindications and Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.

Assorted Dressings, 3 Sizes



AEROPLAST PREMIUM FABRIC ASSORTED DRESSINGS ENV/20

AFP504