

## AEROPLAST PREMIUM FABRIC DRESSING STRIP 7.5CM X 1M

**AFP701** 





#### **ADDITIONAL INFORMATION**

Weight 0.05 kg

**Dimensions**  $3 \times 10.5 \times 9.7 \text{ cm}$ 

Carton Qty 144

Pricing Unit Roll/1

**Size** 7.5cm x 1M

Min Order Qty 12

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311546

**Tax Schedule** Taxable

Item Length 100

Item Width 7.5

Item Measurement Unit cm

**Is Medical Device** Yes

Is Serialized Item

**Is Dangerous Good** No

No

Shelf Life (Months) 60

#### **PRODUCT DESCRIPTION**

The AeroPlast™ Premium Weight Fabric Adhesive Bandage is a top-tier, durable dressing designed for extended wear and maximum protection. Featuring our SupaBond™ extra-strong adhesive, this bandage is the ideal choice when long-lasting, reliable coverage is needed.



### AEROPLAST PREMIUM FABRIC DRESSING STRIP 7.5CM X 1M

**AFP701** 

Aeroplast™ Premium Weight Fabric maximum long wearing durable, fabric for ultimate wound protection

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

wound exudate pad super absorption, reducing skin maceration

Non-stick barrier protects wound with a non-adherent barrier, minimising trauma on removal

#### Unmatched SupaBond™ Adhesive

Engineered exclusively for AeroPlast™, the SupaBond™ adhesive offers unrivaled adhesion that secures firmly across the application site, preventing lift-up or edge curling. SupaBond™ is latex-free and leaves no sticky residue, addressing common issues found with heavy-duty fabric plasters.

**Superior SupaSorb™ Wound Pad** The SupaSorb™ wound pad is designed for optimal wound management, effectively absorbing blood and exudates into its internal pockets. This super-absorbent pad reduces the frequency of dressing changes, protecting the wound and keeping it clean and comfortable. The non-adherent film layer minimizes pain and trauma during dressing changes.

#### **Heavy-Duty Fabric Backing**

Constructed with conforming and air-permeable heavy-duty fabric, the AEROPLAST™ Premium Fabric dressing is perfect for demanding applications, such as in the manufacturing industry. Its superior adhesion and flexibility ensure the bandage stays in place and conforms to body contours, maintaining wound cleanliness and comfort under the toughest conditions.

The 7.5cm x 1m premium fabric strip bandage is designed for medium to large wounds, offering a size that can be cut to size. Offering a durable protective layer to create a moist healing environment.

Available as a roll of 1.

#### **Contraindications and Precautions**

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

**Dressing Strip** 



# AEROPLAST PREMIUM FABRIC DRESSING STRIP 7.5CM X 1M

**AFP701**