



## ADDITIONAL INFORMATION

|                             |                 |
|-----------------------------|-----------------|
| <b>Weight</b>               | 0.022 kg        |
| <b>Dimensions</b>           | 0.2 × 9 × 11 cm |
| <b>Carton Qty</b>           | 375             |
| <b>Pricing Unit</b>         | Env/12          |
| <b>Size</b>                 | Assorted        |
| <b>Min Order Qty</b>        | 15              |
| <b>ARTG ID</b>              | 159198          |
| <b>GMDN Code</b>            | 44990           |
| <b>Medical Device Class</b> | Class 1s        |
| <b>UNSPSC Code</b>          | 42311546        |
| <b>Tax Schedule</b>         | Taxable         |
| <b>Is Medical Device</b>    | Yes             |
| <b>Is Serialized Item</b>   | No              |
| <b>Is Dangerous Good</b>    | No              |
| <b>Shelf Life (Months)</b>  | 60              |



## PRODUCT DESCRIPTION

The AeroPlast™ Premium Weight Fabric Adhesive Bandage is a top-tier, durable dressing designed for extended wear and maximum protection. Featuring our SupaBond™ extra-strong adhesive, this bandage is the ideal choice when long-lasting, reliable coverage is needed.

Aeroplast™  
Premium Weight Fabric  
maximum long wearing durable,  
fabric for ultimate wound protection

Supabond™  
Hypo-allergenic adhesive  
exclusively engineered,  
for low irritation and high adhesion

wound  
exudate pad  
super absorption,  
reducing skin maceration

Non-stick barrier  
protects wound with a non-adherent barrier,  
minimising trauma on removal

### **Unmatched SupaBond™ Adhesive**

Engineered exclusively for AeroPlast™, the SupaBond™ adhesive offers unrivaled adhesion that secures firmly across the application site, preventing lift-up or edge curling. SupaBond™ is latex-free and leaves no sticky residue, addressing common issues found with heavy-duty fabric plasters.

**Superior SupaSorb™ Wound Pad** The SupaSorb™ wound pad is designed for optimal wound management, effectively absorbing blood and exudates into its internal pockets. This super-absorbent pad reduces the frequency of dressing changes, protecting the wound and keeping it clean and comfortable. The non-adherent film layer minimizes pain and trauma during dressing changes.

### **Heavy-Duty Fabric Backing**

Constructed with conforming and air-permeable heavy-duty fabric, the AEROPLAST™ Premium Fabric dressing is perfect for demanding applications, such as in the manufacturing industry. Its superior adhesion and flexibility ensure the bandage stays in place and conforms to body contours, maintaining wound cleanliness and comfort under the toughest conditions.

**Fingertip & Knuckle Dressings:** contoured shapes conform to joints for a secure, flexible fit, supporting strong adhesion on high-motion areas while maintaining breathable, skin-friendly coverage.

Supplied in an **envelope of 12** dressings to streamline kit top-ups and personal carry. This pack helps reduce restocking frequency while ensuring strong adhesion, breathable comfort, and a conforming fit on hands and joints.

### **Contraindications and Precautions**

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



# AEROPLAST PREMIUM FABRIC FINGERTIP & KNUCKLE DRESSINGS ENV/12

AFP312

F/tip & Knuckle Dressings x 12