

AEROPLAST PREMIUM FABRIC FINGERTIP & KNUCKLE DRESSINGS ENV/12

AFP312





ADDITIONAL INFORMATION

Weight 0.022 kg

Dimensions $0.2 \times 9 \times 11$ cm

Carton Qty 375

Pricing Unit Env/12

Size Assorted

Min Order Qty 15

ARTG ID 159198

GMDN Code 44990

Medical Device Class Class 1s

UNSPSC 42311546

Taxable

Tax Schedule Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

The AeroPlast™ Premium Weight Fabric Adhesive Bandage is a top-tier, durable dressing designed for extended wear and maximum protection. Featuring our SupaBond™ extra-strong adhesive, this bandage is the ideal choice when long-lasting, reliable coverage is needed.

Aeroplast™ Premium Weight Fabric maximum long wearing durable,



AEROPLAST PREMIUM FABRIC FINGERTIP & KNUCKLE DRESSINGS ENV/12

AFP312

fabric for ultimate wound protection

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

wound exudate pad super absorption, reducing skin maceration

Non-stick barrier protects wound with a non-adherent barrier, minimising trauma on removal

Unmatched SupaBond™ Adhesive

Engineered exclusively for AeroPlast™, the SupaBond™ adhesive offers unrivaled adhesion that secures firmly across the application site, preventing lift-up or edge curling. SupaBond™ is latex-free and leaves no sticky residue, addressing common issues found with heavy-duty fabric plasters.

Superior SupaSorb™ Wound Pad The SupaSorb™ wound pad is designed for optimal wound management, effectively absorbing blood and exudates into its internal pockets. This super-absorbent pad reduces the frequency of dressing changes, protecting the wound and keeping it clean and comfortable. The non-adherent film layer minimizes pain and trauma during dressing changes.

Heavy-Duty Fabric Backing

Constructed with conforming and air-permeable heavy-duty fabric, the AEROPLAST™ Premium Fabric dressing is perfect for demanding applications, such as in the manufacturing industry. Its superior adhesion and flexibility ensure the bandage stays in place and conforms to body contours, maintaining wound cleanliness and comfort under the toughest conditions.

The AeroPlast™ Premium Weight Fabric Knuckle and Fingertip Dressings are expertly crafted to provide superior protection and comfort for knuckle and fingertip injuries.

An envelope of 12 AeroPlast™ Premium Weight Fabric Adhesive Bandages offers a compact and lightweight option, perfect for on-thego use in personal first aid kits, travel packs, or car glove compartments, ensuring you're prepared with essential dressings without the bulkiness of a box.

Contraindications and Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

F/tip & Knuckle Dressings x 12



AEROPLAST PREMIUM FABRIC FINGERTIP & KNUCKLE DRESSINGS ENV/12

AFP312