

## AEROPLAST PREMIUM FABRIC STANDARD STRIP 7.2 X 1.9CM BOX/50

AFP5015





### **ADDITIONAL INFORMATION**

Weight 0.06 kg

**Dimensions**  $3.3 \times 8.8 \times 10 \text{ cm}$ 

Carton Qty 72

Pricing Unit Box/50

**Size** 7.2cm x 1.9cm

Min Order Qty 12

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311546

**Tax Schedule** Taxable

Item Length 7.2

Item Width 1.9

Item Measurement Unit cm

**Is Medical Device** Yes

**Is Serialized Item** No

**Is Dangerous Good** No

Shelf Life (Months) 60

### **PRODUCT DESCRIPTION**

AEROPLAST Premium Fabric adhesive bandages deliver durable, long-wearing coverage for minor cuts, grazes and abrasions in workplaces and everyday first aid. A premium-weight woven backing conforms to fingers, knuckles and joints while allowing airflow to the skin for comfort. SupaBond high-adhesion technology provides a strong, secure hold that resists edge lift during frequent hand use yet peels away cleanly when removed correctly. The adhesive system is latex-free and hypo-allergenic for broad skin compatibility, supporting use on sensitive skin.

The dressing core features a SupaSorb absorbent pad with a non-adherent contact layer to manage light exudate and minimise trauma at change. Individually wrapped, sterile dressings are cold-sealed with clear pull tabs for fumble-free opening and tamper protection. The breathable, open-weave fabric backing is air-permeable to help reduce moisture build-up under the dressing, while the conforming fit maintains coverage across contoured areas and during movement.



# AEROPLAST PREMIUM FABRIC STANDARD STRIP 7.2 X 1.9CM BOX/50

AFP5015

Built for reliability in busy environments, the heavy-duty fabric resists abrasion and maintains adhesion in challenging conditions. Sterile single-use dressings are ethylene oxide sterilised and clearly labelled with storage guidance (0°C–30°C) and change-daily instructions to support correct use in first aid rooms, cabinets and mobile kits. Efficient dispensing formats help streamline replenishment for facilities and service providers.

- **Strong adhesion:** SupaBond technology holds fast and resists edge curl.
- **Skin compatible:** hypo-allergenic, latex-free adhesive suitable for sensitive skin.
- **Breathable backing:** air-permeable fabric supports comfort during wear.
- **Absorbent, non-stick pad:** manages light exudate and reduces trauma on removal.
- Sterile & easy-open: cold-sealed wrappers with pull tabs for hygienic application.
- Conforming fit: flexible fabric maintains coverage over digits and joints.
- Work-ready durability: premium fabric designed for demanding conditions.

The **Size – Standard Strip 7.2 × 1.9 cm** variant provides everyday coverage for fingers and small cuts, offering conforming fabric, strong adhesion and breathable comfort for extended wear.

**Box of 50** provides efficient replenishment for first aid kits and small workplaces, keeping a practical number on hand while the compact carton and tear-off dispenser support quick, hygienic access.

#### Contraindications and Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Fabric Standard Strip x 50