



## ADDITIONAL INFORMATION

<b>Weight</b>	0.12 kg
<b>Dimensions</b>	4.1 × 12.8 × 10 cm
<b>Carton Qty</b>	36
<b>Pricing Unit</b>	Box/100
<b>Size</b>	7.5cm x 2.5cm
<b>Min Order Qty</b>	12
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC</b>	42311546
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	7.5
<b>Item Width</b>	2.5
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Barcode (GTIN/EAN/UPC)</b>	9341394000537

## PRODUCT DESCRIPTION

The AeroPlast Premium Visual plaster is a premium latex free, hypo-allergenic, visual detectable plaster specifically designed for the retail food industry and non-commercial food preparation



# AEROPLAST PREMIUM VISUAL STRIP 7.5 X 2.5CM BOX/100

AV2000

Advanced  
Film Technology  
Unique Vibrant Blue with Superior  
Stretch and Adhesion

Supabond™  
Hypo-allergenic adhesive  
exclusively engineered,  
for low irritation and high adhesion

SuperSorb™ Wound  
Exudate Pad  
super absorption,  
reducing skin maceration

Non-stick Barrier  
Protects Wound with a Non-Adherent Barrier,  
Minimising Trauma on Removal

## Visually Detectable

- AeroPlast Premium Visual comes in a vibrant blue colour, which is the only colour that does not occur naturally in food, making it easily detectable.

## SupaSorb™ Wound Pad

- An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

## SupaBond™ Adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling. The unique specifications of film and adhesive are key to the products reliability, ensuring that no matter how wet, steam-filled or hot the environment becomes, this product cannot become separated from the wearer.

## Premium plaster backing film

- AeroPlast Visual is made with the ultimate in film technology, providing a bacterial barrier which resists the ingress of water, oil and grease. The standard plaster size is larger for even greater wound isolation.

## Contraindications / Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

**This product is not a substitute for direct medical attention from a healthcare practitioner**

Extra Wide Strip



**AEROPLAST PREMIUM VISUAL STRIP 7.5 X 2.5CM  
BOX/100**

AV2000