

AV2000





ADDITIONAL INFORMATION

0.12 kg
4.1 × 12.8 × 10 cm
36
Box/100
7.5cm x 2.5cm
12
159198
44990
Class 1s
42311546
Taxable
7.5
2.5
cm
Yes
No
No
36

PRODUCT DESCRIPTION

The AeroPlast Premium Visual plaster is a premium latex free, hypo-allergenic, visual detectable plaster specifically designed for the retail food industry and non-commercial food preparation



AV2000

Advanced Film Technology Unique Vibrant Blue with Superior Stretch and Adhesion

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

SuperSorb[™] Wound Exudate Pad super absorption, reducing skin maceration

Non-stick Barrier Protects Wound with a Non-Adherent Barrier, Minimising Trauma on Removal



AV2000

Visually Detectable

• AeroPlast Premium Visual comes in a vibrant blue colour, which is the only colour that does not occur naturally in food, making it easily detectable.

SupaSorb[™] Wound Pad

 An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb[™] leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

SupaBond[™] Adhesive

• SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling. The unique specifications of film and adhesive are key to the products reliability, ensuring that no matter how wet, steam-filled or hot the environment becomes, this product cannot become separated from the wearer.

Premium plaster backing film

• Aeroplast Visual is made with the ultimate in film technology, providing a bacterial barrier which resists the ingress of water, oil and grease. The standard plaster size is larger for even greater wound isolation.

Contraindications / Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Extra Wide Strip



AV2000