



**ADDITIONAL INFORMATION**

<b>Weight</b>	0.12 kg
<b>Dimensions</b>	4.1 × 12.8 × 10 cm
<b>Carton Qty</b>	36
<b>Pricing Unit</b>	Box/100
<b>Size</b>	7.5cm x 2.5cm
<b>Min Order Qty</b>	12
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311546
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	7.5
<b>Item Width</b>	2.5
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	36



**PRODUCT DESCRIPTION**

AEROPLAST Premium Visual is a blue, visually detectable adhesive bandage range engineered for food preparation, catering, and processing environments where fast identification matters. The vibrant blue backing contrasts clearly with food and stainless-steel surfaces to support routine hygiene checks. Each dressing is latex-free and hypo-allergenic to suit diverse teams and sensitive skin while maintaining dependable performance throughout demanding shifts. Sterile, single-use, individually wrapped dressings protect

product integrity until application.

The premium waterproof film forms a barrier to water, oils, and grease while remaining air-permeable for comfort. A conforming, flexible construction bends and stretches with hands and joints to maintain coverage during frequent handwashing and repetitive tasks. Adhesion strength is delivered by SupaBond technology, designed for strong initial tack and long wear. It helps resist edge lift and curling in hot, steamy, or wet conditions common to food service and production lines, supporting continuous protection.

Beneath the film, the SupaSorb embossed wound pad manages exudate by drawing fluid into micro-pockets, helping keep the site cleaner and more comfortable. A non-adherent contact layer minimises trauma on removal to support skin condition during routine changes. The overall design balances robust adhesion, breathability, conformability, and skin compatibility—key attributes for a reliable visual plaster used in professional kitchens and home food preparation.

From workplace first aid cabinets to travel kits, this range provides durable, air-permeable coverage with strong adhesion and skin-friendly removal. It delivers consistent performance that supports hygiene protocols without compromising wearer comfort.

- **Blue visual design:** easy to spot in food areas to aid hygiene checks.
- **Waterproof yet breathable:** barrier to water, oils, and grease with air-permeable comfort.
- **Strong adhesion:** SupaBond technology helps prevent edge lift and curling.
- **Comfort wound pad:** SupaSorb pad with non-adherent layer for gentler removal.
- **Latex-free and hypo-allergenic:** designed for broad skin compatibility.

The **Strip 7.5 cm x 2.5 cm** size conforms to fingers, hands, and wrists. It provides a balance of strong adhesion and flexible comfort, making it well suited to frequent hand movement in food handling, catering, and everyday first aid.

Supplied in a **box of 100** blue visual bandages to streamline bulk replenishment for food areas, workplace first aid stations, and home or travel kits, helping reduce mid-shift changeovers and keeping high-use locations well stocked.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Extra Wide Strip