

### **AEROPLAST PREMIUM VISUAL STRIP 7.5 X 2.5CM BOX/50**

AV2050





### **ADDITIONAL INFORMATION**

Weight 0.065 kg

**Dimensions**  $4.1 \times 8 \times 10.2 \text{ cm}$ 

**Carton Qty** 60

**Pricing Unit** Box/50

Size 7.5cm x 2.5cm

**Min Order Qty** 12

**ARTG ID** 159198 **GMDN Code** 44990

**Medical Device Class** Class 1s

UNSPSC 42311546

**Tax Schedule** Taxable

7.5 **Item Length Item Width** 2.5

Item Measurement Unit cm

Is Serialized Item

No

36

**Is Dangerous Good** No

Shelf Life (Months)

#### **PRODUCT DESCRIPTION**

The AeroPlast Premium Visual plaster is a premium latex free, hypo-allergenic, visual detectable plaster specifically designed for the retail food industry and non-commercial food preparation



## AEROPLAST PREMIUM VISUAL STRIP 7.5 X 2.5CM BOX/50

AV2050

Advanced Film Technology Unique Vibrant Blue with Superior Stretch and Adhesion

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

SuperSorb™ Wound Exudate Pad super absorption, reducing skin maceration

Non-stick Barrier Protects Wound with a Non-Adherent Barrier, Minimising Trauma on Removal



# AEROPLAST PREMIUM VISUAL STRIP 7.5 X 2.5CM BOX/50

AV2050

### **Visually Detectable**

• AeroPlast Premium Visual comes in a vibrant blue colour, which is the only colour that does not occur naturally in food, making it easily detectable.

Extra Wide Strip