

### AEROPLAST PREMIUM VISUAL STRIP 7.5 X 2.5CM BOX/50

AV2050





#### **ADDITIONAL INFORMATION**

Weight 0.065 kg

**Dimensions**  $4.1 \times 8 \times 10.2 \text{ cm}$ 

Carton Qty 60

Pricing Unit Box/50

**Size** 7.5cm x 2.5cm

Min Order Qty 12

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311546

**Tax Schedule** Taxable

Item Length7.5Item Width2.5

Item Measurement Unit cm

Is Medical Device Yes

**Is Serialized Item** No

**Is Dangerous Good** No

Shelf Life (Months) 36

#### **PRODUCT DESCRIPTION**

The AeroPlast Premium Visual plaster is a premium latex free, hypo-allergenic, visual detectable plaster specifically designed for the retail food industry and non-commercial food preparation



# AEROPLAST PREMIUM VISUAL STRIP 7.5 X 2.5CM BOX/50

AV2050

Advanced Film Technology Unique Vibrant Blue with Superior Stretch and Adhesion

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

SuperSorb™ Wound Exudate Pad super absorption, reducing skin maceration

Non-stick Barrier Protects Wound with a Non-Adherent Barrier, Minimising Trauma on Removal



# AEROPLAST PREMIUM VISUAL STRIP 7.5 X 2.5CM BOX/50

AV2050

### **Visually Detectable**

• AeroPlast Premium Visual comes in a vibrant blue colour, which is the only colour that does not occur naturally in food, making it easily detectable.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Extra Wide Strip