



### ADDITIONAL INFORMATION

<b>Weight</b>	0.05 kg
<b>Dimensions</b>	3.2 × 8.5 × 10 cm
<b>Carton Qty</b>	120
<b>Pricing Unit</b>	Box/20
<b>Size</b>	Assorted
<b>Min Order Qty</b>	6
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311546
<b>Tax Schedule</b>	Taxable
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60



### PRODUCT DESCRIPTION

AeroPlast Sensitive provides gentle, effective wound coverage for delicate, allergy-prone, or compromised skin. Each dressing combines a soft, breathable spunlace non-woven backing with a low-trauma contact layer and an absorbent wound pad for dependable daily protection. Dressings are individually wrapped in a paper wrapper to protect sterility and product integrity until use.

The SupaBond adhesive offers a strong yet skin-kind hold that resists edge lift while supporting comfortable wear on hands, joints, and high-movement areas. The formulation is hypo-allergenic and latex-free, helping reduce the likelihood of irritation on sensitive or fragile skin.

The non-adherent barrier helps prevent the pad from sticking to the wound, supporting low-pain removal. The air-permeable fabric

backing maintains breathability and conforms smoothly to curves and contours for reliable coverage. An absorbent, non-woven sterile wound pad manages light exudate to help keep the area cleaner between changes.

- **Gentle, secure hold:** SupaBond adhesive balances adhesion strength with skin comfort.
- **Skin-friendly design:** hypo-allergenic and latex-free for sensitive or fragile skin.
- **Air-permeable backing:** conforming non-woven fabric supports comfort during wear.
- **Low-trauma removal:** non-adherent contact layer reduces sticking to the wound.
- **Individually wrapped & sterile:** EO-sterilised plasters protected in paper wrappers.

Assorted sizes in a box of **20** provide versatile coverage for fingers, hands, and small areas, streamlining replenishment and supporting quick selection for everyday minor cuts and grazes.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.